

Bridgeport Pathology IHC

GEN NAME	CLONE	DILUTION	AR	MANUFACTURER	CATALOG #	PLATFORM
AE1/AE3 (CK)	AE1/AE3	Predilute	TRS Lo. 20'	Dako	IR05361-2	Link 48
AE1/AE3 (CK)	AE1/AE3	Predilute	TRS Hi 30'	Dako	GA05361-2	Omnis
AFP	Rab Poly	Predilute	TRS Hi 20'	Dako	IR50061-2	Link 48
ANDR REC	AR 441	1:50	TRS Hi 20'	Dako	MAB8052	Link 48
ANDR REC	AR 442	1:25	ER2 20'	Dako	PA0963	Bond
ANDR REC	AR 442	1:50	TRS Hi 30'	Dako	AFP-L-CE	Omnis
Arginase	EP261	Predilute	TRS Hi 20'	BioCare	API3058AA	Link 48
Arginase	EP261	Predilute	TRS Hi 30'	BioCare	API3058AA	Omnis
BAP1	C-4	1:50	TRS Hi 20'	Santa Cruz	sc-28383	Link 48
BAP1	C-4	1:50	TRS Hi 30'	Santa Cruz	sc-28383	Omnis
BCL-2	123	Predilute	TRS Hi 20'	Dako	IR61461-2	Link 48
BCL-6	PG-B6p	Predilute	TRS Lo. 30'	Dako	GA62561-2	Omnis
Beta Cat.	β-Catenin-1	Predilute	TRS Hi 20'	Dako	IR70261-2	Link 48
Beta Cat.	β-Catenin-1	Predilute	TRS Lo 30'	Dako	GA70261-2	Omnis
Calcitonin	Rab Poly	Predilute	TRS Hi 20'	Dako	IR51561-2	Link 48
Calcitonin	Rab Poly	Predilute	TRS Lo 30'	Dako	GA51561-2	Omnis
Caldesmon	h-CD	Predilute	TRS Hi 20'	Dako	IR05461-2	Link 48
Calponin	CALP	1:100	TRS Hi 20'	Dako	M355601-2	Link 48
Calponin	CALP	Predilute	TRS Hi 30'	BioCare	PM172AA	Omnis
Calretinin	DAK CALRET 1	Predilute	TRS Hi 20'	Dako	IR62761-2	Link 48
Calretinin	DAK CALRET 1	Predilute	TRS Hi 30'	Dako	IR62761-2	Omnis
CAM5.2	CAM5.2	1:5	TRS Lo 20'	BD	349205	Link 48
CAIX	EP161	Predilute	TRS Hi 30'	Cell Marque	379R-18	Omnis
CD1a	O10	Predilute	TRS Lo 30'	Dako	IR06961-2	Omnis
CD2	AB75	Predilute	TRS Hi 30'	Dako	GA65161-2	Omnis
CD3	Rab Poly	Predilute	TRS Hi 20'	Dako	IR50361-2	Link 48
CD3	Rab Poly	Predilute	TRS Lo 30'	Dako	GA50361-2	Omnis
CD4	4B12	Predilute	TRS Hi 30'	Dako	IR64961-2	Omnis
CD5	4C7	Predilute	TRS Hi 20'	Dako	IR08261-2	Link 48
CD5	4C7	Predilute	TRS Hi 30'	Dako	IR08261-2	Omnis
CD8	C8/144B	Predilute	TRS Hi 30'	Dako	GA62362-1	Omnis
CD10	56C6	Predilute	TRS Hi 20'	Dako	IR64861-2	Link 48
CD10	56C6	Predilute	TRS Hi 30'	Dako	GA64861-2	Omnis
CD15	Carb-3	Predilute	TRS Hi 20'	Dako	IR06261-2	Link 48
CD19	LE-CD19	Predilute	TRS Hi 30'	Dako	GA65661-2	Omnis
CD20	L26	Predilute	TRS Hi 20'	Dako	IR60461-2	Link 48
CD20	L26	Predilute	TRS Hi 30'	Dako	GA60461-2	Omnis
CD21	1F8	Predilute	TRS Hi 30'	Dako	IR60861-2	Omnis
CD23	DAK-CD23	Predilute	TRS Hi 20'	Dako	IR78161-2	Link 48
CD23	DAK-CD23	Predilute	TRS Hi 30'	Dako	GA78161-2	Omnis

GEN NAME	CLONE	DILUTION	AR	MANUFACTURER	CATALOG #	PLATFORM
CD30	BER-H2	Predilute	TRS Hi 20'	Dako	IR60261-2	Link 48
CD30	BER-H2	Predilute	TRS Hi 30'	Dako	GA60162-2	Omnis
CD31	JC70A	Predilute	TRS Hi 20'	Dako	IR61061-2	Link 48
CD31	JC70A	Predilute	TRS Hi 30'	Dako	GA61061-2	Omnis
CD34	QBEnd-10	Predilute	TRS Hi 20'	Dako	IR63261-2	Link 48
CD34	QBEnd-10	Predilute	TRS Hi 30'	Dako	GA63261-2	Omnis
CD38	SP 149	1:100	TRS Hi 20'	Cell Marque	118R-16	Link 48
CD38	SP 149	1:100	ER1 20'	Cell Marque	118R-16	Bond
CD38	SP 149	1:100	TRS Hi 30'	Cell Marque	118R-16	Omnis
CD43	DF-T1	Predilute	TRS Hi 20'	Dako	IR63661-2	Link 48
CD45 (LCA)	2B11+PD7/26	Predilute	TRS Hi 20'	Dako	IR75161-2	Link 48
CD45 (LCA)	2B11+PD7/26	Predilute	TRS Hi 30'	Dako	GA75161-2	Omnis
CD56	123C3	Predilute	TRS Hi 20'	Dako	IR62861-2	Link 48
CD56	123C3	Predilute	TRS Hi 30'	Dako	IR62861-2	Omnis
CD68 (KPI)	KP1	Predilute	TRS Hi 20'	Dako	IR60961-2	Link 48
CD68 (PGM1)	PG-M1	Predilute	TRS Hi 30'	Dako	GA61361-2	Omnis
CD79a	JCB117	Predilute	TRS Hi 30'	Dako	GA62161-2	Omnis
CD99	12E7	Predilute	TRS Hi 20'	Dako	IR05761-2	Link 48
CD99	12E7	Predilute	TRS Hi 30'	Dako	IR05761-2	Omnis
CD117 (C-Kit)	YR145	1:400	TRS Hi 20'	Cell Marque	117R-16-ASP	Link 48
CD117 (C-Kit)	EP10	Predilute	TRS Hi 30'	BioCare	PME296AA	Omnis
CD138	MI15	Predilute	TRS Hi 20'	Dako	IR64261-2	Link 48
CD138	B-I38	Predilute	TRS Hi 30'	Dako	GA64261-2	Omnis
CDX2	DAK-CDX2	Predilute	TRS Hi 20'	Dako	IR08061-2	Link 48
CDX2	DAK-CDX2	Predilute	TRS Hi 30'	Dako	GA08061-2	Omnis
CDX2/CK20	AK-CDX2/Ks20.8	Predilute	TRS Hi 20'	Dako	IR61-2/IR77761-2	Link 48
CDX2/CK20	AK-CDX2/Ks20.8	Predilute	TRS Hi 30'	Dako	GA71-2/GA77761-2	Omnis
CEA M	II-7	Predilute	TRS Hi 20'	Dako	IR62261-2	Link 48
CEA M	II-7	Predilute	None	Dako	GA62261-2	Bond
CEA M	II-7	Predilute	TRS Hi 30'	Dako	GA62261-2	Omnis
CEA P	Rab Poly	Predilute	TRS Hi 20'	Cell Marque	236A-18	Link 48
CEA P	Rab Poly	Predilute	ER2 5'	Cell Marque	236A-18	Bond
CK5/6	D5/16 B4	Predilute	TRS Hi 20'	Dako	IR78061-2	Link 48
CK5/6	D5/16 B4	Predilute	TRS Hi 30'	Dako	GA78061-2	Omnis
CK7	OV.TL 12/30	Predilute	TRS Hi 20'	Dako	IR61961-2	Link 48
CK7	OV.TL 12/30	Predilute	TRS Hi 30'	Dako	GA61961-2	Omnis
CK19	RCK108	Predilute	TRS Hi 20'	Dako	IR61561-2	Link 48
CK19	RCK108	Predilute	TRS Hi 30'	Dako	GA61561-2	Omnis
CK20	Ks20.8	Predilute	TRS Hi 20'	Dako	IR77761-2	Link 48
CK20	Ks20.8	Predilute	TRS Hi 30'	Dako	GA77761-2	Omnis
Ker903 (CKHMW)	34BE12	Predilute	TRS Hi 30'	Dako	GA05161-2	Omnis
Chromo-A	LK2H10+PHE5	Predilute	TRS Hi 20'	BioCare	PM010AA	Link 48

GEN NAME	CLONE	DILUTION	AR	MANUFACTURER	CATALOG #	PLATFORM
Chromo-A	LK2H10+PHE5	Predilute	TRS Hi 30'	BioCare	PM010AA	Omnis
CMV	DT10+BC90	Predilute	TRS Hi 20'	BioCare	APA118AA	Link 48
CMV	DT10+BC90	Predilute	TRS Hi 30'	BioCare	APA118AA	Omnis
Cyclin-D1	EP12	Predilute	TRS Hi 30'	Dako	GA08361-2	Omnis
D2-40 (Podoplan	D2-40	Predilute	TRS Hi 20'	Dako	IR07261-2	Link 48
D2-40 (Podoplan	D2-40	Predilute	TRS Hi 30'	Dako	IR07261-2	Omnis
Desmin	D33	Predilute	TRS Hi 20'	Dako	IR60661-2	Link 48
Desmin	D33	Predilute	TRS Hi 30'	Dako	IR60661-2	Omnis
E-CAD	NCH-38	Predilute	TRS Hi 20'	Dako	IR05961-2	Link 48
E-CAD	NCH-38	Predilute	TRS Hi 30'	Dako	GA05961-2	Omnis
EMA	E29	Predilute	TRS Hi 20'	Dako	IR62961-2	Link 48
EMA	E29	Predilute	TRS Hi 30'	Dako	GA62961-2	Omnis
ER (SP1)	SP1	1:100	TRS Hi 20'	Epedia	RM-9101-S	Link 48
ERG	EP111	Predilute	TRS Hi 20'	Dako	IR65961-2	Link 48
ERG	EP111	Predilute	TRS Hi 30'	Dako	GA65961-2	Omnis
Gastrin	Rab Poly	Predilute	TRS Hi 30'	Dako	GA51961-2	Omnis
Gastrin	Rab Poly	Predilute	TRS Hi 20'	Dako	IR51961-2	Link 48
GATA3	L50-823	Predilute	TRS Hi 20'	BioCare	PM405AA	Link 48
GATA3	L50-823	Predilute	TRS Hi 30'	BioCare	PM405AA	Omnis
GFAP	Rab Poly	Predilute	TRS Hi 30'	Dako	GA52461-2	Omnis
Glypican	IG12	1:100	TRS Hi 20'	Cell Marque	261A-96	Link 48
Glypican	IG12	Predilute	TRS Hi 30'	BioCare	PM396AA	Omnis
HCG	Rab Poly	1:200	TRS Hi 20'	Cell Marque	234A-16	Link 48
HCG	Rab Poly	1:200	ER1 10'	Cell Marque	234A-16	Bond
HEP-PAR1	OCHIES	Predilute	TRS Hi 20'	Dako	IR62461-2	Link 48
HEP-PAR1	OCHIES	Predilute	TRS Hi 30'	Dako	GA62461-2	Omnis
Her2/Neu	EP3	1:100	TRS Hi 20'	Cell Marque	AC0014	Link 48
HMB-45	HMB-45	Predilute	TRS Hi 20'	Dako	IR05261-2	Link 48
HMB-45	HMB-45	Predilute	TRS Hi 30'	Dako	GA05261-2	Omnis
HP	Rab Poly	1:50	TRS Hi 20'	Dako	B047101-1	Link 48
HP	Rab Poly	1:50	ER2 15'	Dako	B047101-1	Bond
HP	Rab Poly	1:50	TRS Hi 30'	Dako	B047101-1	Omnis
HSV1/II	Rab Poly	Predilute	TRS Hi 20'	BioCare	PP108AA	Link 48
HSV1/II	Rab Poly	Predilute	None	BioCare	PP108AA	Bond
HSV1/II	Rab Poly	Predilute	TRS Hi 30'	BioCare	PP108AA	Omnis
INSM1	A8	1:200	TRS Hi 20'	Santa Cruz	sc-271408	Link 48
INSM1	A8	1:200	TRS Hi 30'	Santa Cruz	sc-271408	Omnis
Inhibin	R1	Predilute	TRS Hi 20'	Dako	IR05861-2	Link 48
Inhibin	R1	Predilute	TRS Hi 30'	Dako	GA05861-2	Omnis
Ki-67	MIB-1	Predilute	TRS Hi 20'	Dako	IR62661-2	Link 48
Ki-67	MIB-1	Predilute	TRS Hi 30'	Dako	GA62661-2	Omnis
Mel-A	A103	Predilute	TRS Hi 20'	Dako	IR63361-2	Link 48

GEN NAME	CLONE	DILUTION	AR	MANUFACTURER	CATALOG #	PLATFORM
Mel-A	A103	Predilute	TRS Hi 30'	Dako	IR63361-2	Omnis
MiTF	34CA5	Predilute	TRS Hi 30'	BioCare	PM423H	Omnis
MLH1	ES05	Predilute	TRS Hi 20'	Dako	IR07961-2	Link 48
MLH1	ES05	Predilute	TRS Hi 30'	Dako	IR07961-2	Omnis
MUSA	HHF-35	Predilute	TRS Hi 20'	Dako	IR70061-2	Link 48
MUSA	HHF-35	Predilute	TRS Hi 30'	Dako	IR70061-2	Omnis
MSH-2	FE11	Predilute	TRS Hi 20'	Dako	IR08561-2	Link 48
MSH-2	FE11	Predilute	TRS Hi 30'	Dako	IR08561-2	Omnis
MSH-6	EP49	Predilute	TRS Hi 20'	Dako	IR08661-2	Link 48
MSH-6	EP49	Predilute	TRS Hi 30'	Dako	IR08661-2	Omnis
MUM1	MUM1p	Predilute	TRS Hi 30'	Dako	GA64461-2	Omnis
Myeloperoxidase	Rab Poly	Predilute	TRS Hi 30'	Dako	GA651161-2	Omnis
Napsin-A	IP64	1:100	TRS Hi 20'	Leica		Link 48
Napsin-A	Rab Poly	Predilute	TRS Hi 30'	BioCare	PP434AA	Omnis
NKX3.1	Rab Poly	Predilute	TRS Hi 20'	BioCare	PP422AA	Link 48
NKX3.1	Rab Poly	Predilute	TRS Hi 30'	BioCare	PP422AA	Omnis
NSE	BBS/NC/VI-H14	Predilute	TRS Hi 30'	Dako	IR61261-2	
PanCK (MNF116)	MNF116	1:800	TRS Hi 20'	Dako	M082101-2	Link 48
PanCK (MNF116)	MNF116	1:800	TRS Hi 30'	Dako	M082101-2	Omnis
PIN4	be12,4A4,13H4	Predilute	TRS Hi 20'	Dako/Biocare	DAA/GA06061-2	Link 48
PIN4	be12,4A4,13H4	Predilute	TRS Hi 30'	Dako/Biocare	DAA/GA06061-2	Omnis
P16	EGH4	1:10	TRS Hi 20'	Ventana	6680003001	Link 48
P40	BC28	1:100	TRS Hi 20'	BioCare	ACI3066C	Link 48
P40	BC28	Predilute	TRS Hi 30'	BioCare	API3066H	Omnis
P40/CK5/6	BC28,D5/16 B4	or Predilute	TRS Hi 20'	Dako/Biocare	66C/GA78061-2	Link 48
P40/CK5/6	BC28,D5/16 B4	or Predilute	TRS Hi 30'	Dako/Biocare	66C/GA78061-2	Omnis
p504s (AMACR)	13H4	Predilute	TRS Hi 30'	Dako	GA06061-2	Omnis
P53	DO-7	Predilute	TRS Hi 20'	Dako	IR61661-2	Link 48
P53	DO-7	Predilute	TRS Hi 30'	Dako	GA61661-2	Omnis
P57	25-B2	1:25	TRS Hi 20'	Leica	NCL-L-p57	Link 48
P57	25-B2	1:25	ER2 15'	Leica	NCL-L-p57	Bond
P57	25-B2	1:25	TRS Hi 30'	Leica	NCL-L-p57	Omnis
Pax-5 (BSAP)	DAK-Pax5	Predilute	TRS Hi 30'	Dako	GA65061-2	Omnis
Pax-8	Rab Poly	Predilute	TRS Hi 20'	BioCare	PP379AA	Link 48
Pax-8	Rab Poly	Predilute	TRS Hi 30'	BioCare	PP379AA	Omnis
PLAP	8A9	Predilute	TRS Hi 20'	Dako	IR77961-2	Link 48
PMS-2	EP51	Predilute	TRS Hi 20'	Dako	IR08761-2	Link 48
PMS-2	EP51	Predilute	TRS Hi 30'	Dako	IR08761-2	Omnis
PR	PgR636	1:50	TRS Hi 20'	Dako	M356901-2	Link 48
PSAP	PASE/4LJ	1:2500	TRS Hi 20'	Dako	M079201-2	Link 48
PSAP	PASE/4LJ	1:2500	TRS Hi 30'	Dako	M079201-2	Omnis
PSA (PSMA)	3E6	Predilute	TRS Hi 20'	Dako	IR08961-2	Link 48



