

Radiology Conferences

Radiology residents will meet with the chief pathology resident:

1. For education about case
2. Be shown how to take images –for submission of case to conference.
3. Fill out form for whole slide scanning request
4. Provide 2-3 key radiology images in jpeg format to correspond with the pathology case.
5. Complete any details needed on Pathology scanning form request

Pathology Resident

1. Will send glass slides and the request forms to Path Digital Imaging via the *slide scanning box* in Surgpath
 - Radiology residents will provide jump drive for file transfer; transfer to iornkeys must be done in person. (estimate images 1GB – 3 GB per image at 40 x magnification)
2. Request De-identified report from Brian Daley, send to radiology resident.
3. Locate gross image in Pathology Image repository, send to radiology resident.

Radiology residents:

- Submit online to conference:
 1. Conference form with history
 2. Gross image .jpeg
 3. Photomicrographs .jpeg
 4. Photomicrographs .jpeg with annotations
- At conference they provide .svs file scanned at 40x

PathDigital

- Scan slides at 40x
- Fill out database information from form
- Place on Spectrum Server
- Send email to radiology resident when scans are complete and files transferred to jump drive.
- Return glass slides to Surg path

Radiology Resident will pick up images at Pathology Digital Imaging and provide corresponding radiology images as jpg files. The chief radiology resident will be contacted if this is not complete.